
PK-53 MOLDED CARRYING CASE



NOTES

- 1. Some of the cut-out inserts may need to be removed or pushed down to t your instrumemts.
- 2. To better protect your instruments, the cut-outs have been designed for tight t.
- 3. The upper lid contains extra storage for your accessories and paper work.
- 4. If you need additional help with your new case, please call us at 1-855-STATIC1 (1-855-782-8421).

The PK-53 has a retractable handle and wheels. The case also contains foam inserts for protection of test equipment and accessories.



© Copyright PROSTAT® Corporation 2013. All rights reserved.

Rev. 2: 4-2014

